

Labycare-Weighing Scales Digital Column Type (SKSS-HS46)



Weighing Scales Digital Column Type Manufacturer by LABYCARE

Description :

Weighing Scales Digital Column Type two in one weighing and height scale

The scale has the advantages of large weighing range, accuracy and high sensitivity.

It can go on working with the battery if the electric power can not be supplied, and the battery is rechargeable with electric power.

Max. weighing: 200kg

Min. value per division: 100g

Height range to be measured: 70-190cm

Min. value of height per division: 0.5cm

Size of packing box: 1000×320×260mm

Load platform area (L×W): 375×275mm

Measurement (L×W×H): 535×275×940mm

We are Manufacturer, Exporters and Suppliers of Weighing Scales Digital Column Type at wholesale competitive prices from India.

Specifications :

Weighing Scales Digital Column Type two in one weighing and height scale

The scale has the advantages of large weighing range, accuracy and high sensitivity.

It can go on working with the battery if the electric power can not be supplied, and the battery is rechargeable with electric power.

Max. weighing: 200kg

Min. value per division: 100g

Height range to be measured: 70-190cm

Min. value of height per division: 0.5cm

Size of packing box: 1000×320×260mm

Load platform area(L×W): 375×275mm

Measurement(L×W×H): 535×275×940mm

We are Manufacturer, Exporters and Suppliers of Weighing Scales Digital Column Type at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com