

Labycare-Pulse Oximeter Handheld (SKSS-IE06)



Pulse Oximeter Handheld Manufacturer by LABYCARE

Description :

Pulse Oximeter is used to measure the oxygen level in blood and heart rate for people with COPD, asthma, Congestive Heart Failure (CHF) .

Bright 7 Segment Display

Durable, Compact and Light weight Pulse Oximeter handheld

Measure Spo2 And Pulse rate

Suitable for Adult, Pediatric and Neonate

With Rechargeable Battery & Charger

TECHNICAL PARAMETERS

Display format: LCD Back light : LCD

SP02 : 35-99%

Pulse ratio : 30-250BPM

Pulse wave : LCD graphic display

Battery voltage low indication : LCD Signals.

Power requirement : Two AAA size 1.5 Alkali battery

We are Exporters and Suppliers of Pulse Oximeter Handheld at wholesale competitive prices from India.

Specifications :

Pulse Oximeter is used to measure the oxygen level in blood and heart rate for people with COPD, asthma, Congestive Heart Failure (CHF) .

Bright 7 Segment Display

Durable, Compact and Light weight Pulse Oximeter handheld

Measure Spo2 And Pulse rate

Suitable for Adult, Pediatric and Neonate

With Rechargeable Battery & Charger

TECHNICAL PARAMETERS

Display format: LCD Back light : LCD

SP02 : 35-99%

Pulse ratio : 30-250BPM

Pulse wave : LCD graphic display

Battery voltage low indication : LCD Signals.

Power requirement : Two AAA size 1.5 Alkali battery

We are Exporters and Suppliers of Pulse Oximeter Handheld at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com