

Labycare-Fetal Monitor 7" Display (SKSS-IE29)



Fetal Monitor 7" Display Manufacturer by LABYCARE

Description :

Fetal Monitor used to continuously record the heartbeat of the fetus and the contractions of the woman's uterus during labor.

Fetal Monitor 7" wide TFT color LCD display

Up to 150 patients data saving

Fetal heart sound play and record in PC

Rotation screen

Water proof Probe

GE / Philips Paper Compatible

Interpretation Function

Multiple language support

Built-in quick guide (LCD Type) Water proof probe

Twins Probe Available

Power Consumption : <20W

Dimension : 96 mm (H) X 326 mm (L) X 276 mm (D)

Net Weight : 5.5 KG

Working Temperature :10-40 C, Relative Humidity 10-80%

External Link : RS-232C

Trend : USB Data Saving

PC Interface (Central Monitoring System)

Optional : Rechargeable Battery

We are Exporters and Suppliers of Fetal Monitor 7" Display with Twin Probe at wholesale competitive prices from India.

Specifications :

Fetal Monitor used to continuously record the heartbeat of the fetus and the contractions of the woman's uterus during labor.

Fetal Monitor 7" wide TFT color LCD display

Up to 150 patients data saving

Fetal heart sound play and record in PC

Rotation screen

Water proof Probe

GE / Philips Paper Compatible

Interpretation Function

Multiple language support

Built-in quick guide (LCD Type) Water proof probe

Twins Probe Available

Power Consumption : <20W

Dimension : 96 mm (H) X 326 mm (L) X 276 mm (D)

Net Weight : 5.5 KG

Working Temperature :10-40 C, Relative Humidity 10-80%

External Link : RS-232C

Trend : USB Data Saving

PC Interface (Central Monitoring System)

Optional : Rechargeable Battery

We are Exporters and Suppliers of Fetal Monitor 7" Display with Twin Probe at wholesale competitive

prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com