

Labycare-Fetal Monitor 8.4" Display (SKSS-IE30)



Fetal Monitor 8.4" Display Manufacturer by LABYCARE

Description :

Fetal Monitor used to continuously record the heartbeat of the fetus and the contractions of the woman's uterus during labor.

Fetal Monitor 8.4"color TFT LCD large Screen

The world's advanced DSP technology Double event marks

Three-in-one probe, reduce the line winding.Code-operation ,easy and convenient

12 Hours Data & Waveform Storage and review

Twins monitoring Function Option.

Built in PC Network can be connected to the central monitoring system. (Optional)

TECHNICAL PARAMETERS

Probe : 9 wafer, wide wave Beam, Pulse working mode high sensitivity

Ultrasound Intensity : <5mW / cm²

Recorder : Thermal Printer

Recorder Paper : Z-fold, thermal paper

Printing width : 112m.m.

Toco Range : 0-100 (relative unit)

Display : FHR, TOCO, FM, Mark ,Time ,Freeze, Replay

Power : AC 220V, 50HZ

Power Consumption : <20W

Dimension : 320mm (L) X 260mm (W) X 80mm (H)

Net Weight : 3.2 KG

We are Exporters and Suppliers of Fetal Monitor 8.4" Display at wholesale competitive prices from India.

Specifications :

Fetal Monitor used to continuously record the heartbeat of the fetus and the contractions of the woman's uterus during labor.

Fetal Monitor 8.4"color TFT LCD large Screen

The world's advanced DSP technology Double event marks

Three-in-one probe, reduce the line winding.Code-operation ,easy and convenient

12 Hours Data & Waveform Storage and review

Twins monitoring Function Option.

Built in PC Network can be connected to the central monitoring system. (Optional)

TECHNICAL PARAMETERS

Probe : 9 wafer, wide wave Beam, Pulse working mode high sensitivity

Ultrasound Intensity : <5mW / cm²

Recorder : Thermal Printer

Recorder Paper : Z-fold, thermal paper

Printing width : 112m.m.

Toco Range : 0-100 (relative unit)

Display : FHR, TOCO, FM, Mark ,Time ,Freeze, Replay

Power : AC 220V, 50HZ

Power Consumption : <20W

Dimension : 320mm (L) X 260mm (W) X 80mm (H)

Net Weight : 3.2 KG

We are Exporters and Suppliers of Fetal Monitor 8.4" Display at wholesale competitive prices from

India.

labycare

www.labycare.com | Email: labycare@gmail.com