

Labycare-Cautery Machine 300 Watt Digital (SKSS-IE43)



Cautery Machine 300 Watt Digital Manufacturer by LABYCARE

Description :

Cautery Machine is also known as Electrocautery Machine and Electro Surgical unit

Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

Cautery Machine 300 Watt Digital is micro controller based model with 6 programs

Power Digital : 300 W

Programs : 6

Capability : All surgeries but till 300 W power

Ideal where 400 W not required

Max Power : 300W

Display : 7 Segment Digital

Programs : 6

Monopolar Cut Modes : Pure Cut, Blend1, Blend2, Endo Cut

Monopolar Coagulation Modes : Spray, Dessicate, Fulgurate

Bipolar Modes : Micro, Macro

Tissue Response Circuit : Present

Patient Plate Monitoring System : Present

Dimension (WxHxD) & Weight : 4.5" x 10" x 11" & 3.8Kg.

Specifications :

Cautery Machine is also known as Electrocautery Machine and Electro Surgical unit

Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

Cautery Machine 300 Watt Digital is micro controller based model with 6 programs

Power Digital : 300 W

Programs : 6

Capability : All surgeries but till 300 W power

Ideal where 400 W not required

Max Power : 300W

Display : 7 Segment Digital

Programs : 6

Monopolar Cut Modes : Pure Cut, BlendI, Blend2, Endo Cut

Monopolar Coagulation Modes : Spray, Dessicate, Fulgurate

Bipolar Modes : Micro, Macro

Tissue Response Circuit : Present

Patient Plate Monitoring System : Present

Dimension (WxHxD) & Weight : 4.5" x 10" x 11" & 3.8Kg.

labycare

www.labycare.com | Email: labycare@gmail.com