

Labycare-Cautery Machine 400 Watt Analog (SKSS-IE)



Cautery Machine 400 Watt Analog Manufacturer by LABYCARE

Description :

Cautery Machine is also known as Electrocautery Machine

Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

Cautery Machine is a 400W analog model with number of cut, blend and coagulation modes to do under water surgery.

Power : 400 W Analog Non Programmable

Capability : All surgeries

Robust quality : Great price

Max Power : 400W

Programs : NA

Monopolar Cut Modes : Pure Cut, Blend1, Blend2

Monopolar Coagulation Modes : Spray, Dessicate, Fulgurate

Bipolar Modes : Coag

Tissue Response Circuit : NA

Patient Plate Monitoring System : Present

Dimension (WxHxD) & Weight : 4.5" x 10" x 11" & 3.8 Kg.

Specifications :

Cautery Machine is also known as Electrocautery Machine

Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

Cautery Machine is a 400W analog model with number of cut, blend and coagulation modes to do under water surgery.

Power : 400 W Analog Non Programmable

Capability : All surgeries

Robust quality : Great price

Max Power : 400W

Programs : NA

Monopolar Cut Modes : Pure Cut, Blend1, Blend2

Monopolar Coagulation Modes : Spray, Dessicate, Fulgurate

Bipolar Modes : Coag

Tissue Response Circuit : NA

Patient Plate Monitoring System : Present

Dimension (WxHxD) & Weight : 4.5" x 10" x 11" & 3.8 Kg.

labycare

www.labycare.com | Email: labycare@gmail.com