

Labycare-halogen light source (SKSS-IE49)



halogen light source

Description :

Halogen Light Source is used for High Intensity Halogen Light

24V, 250W excellent distortion bulb

Twin lamp outlet with easy switching

Continuous dimming from 0 to 100%

Easy to replace bulb even if stacked (fused)

Heavy duty transformer ensures uniform and standard power supply

And avoids frequent bulb failure

Fuse blown Indicator with Spare fuse

Customer can replace lamp and fuse easily

We are Exporters and Suppliers of Halogen Light Source at wholesale competitive prices from India.

Specifications :

Halogen Light Source is used for High Intensity Halogen Light

24V, 250W excellent distortion bulb

Twin lamp outlet with easy switching

Continuous dimming from 0 to 100%

Easy to replace bulb even if stacked (fused)

Heavy duty transformer ensures uniform and standard power supply

And avoids frequent bulb failure

Fuse blown Indicator with Spare fuse

Customer can replace lamp and fuse easily

We are Exporters and Suppliers of Halogen Light Source at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com