

Labycare-Vessel Sealer (SKSS-IE51)



Vessel Sealer Manufacturer by LABYCARE

Description :

Electrothermal Bipolar Vessel Sealing (EBVS) is an electrosurgery technology for sealing blood vessels of up to 7mm in diameter

Electro Skin Cautery is used in Electrosurgery

Electrosurgery is done to destroy benign and malignant lesions, to control bleeding, and to cut or excise tissue.

Standard Vessel Sealing Cum Electrosurgical unit.

Graphical LCD Display.

400w Power.

99 programmable.

Patient Plate Conductivity monitoring system.

99 Programs Memory

Display : Graphical LCD Display

Dimension (W x H x D) : 360mm x 145 mm x 305 mm

Weight : 4.5 kg

Input Supply Voltage : 230 VAC + 15%@ 50/60Hz (110 V optional)

Specifications :

Electrothermal Bipolar Vessel Sealing (EBVS) is an electrosurgery technology for sealing blood vessels of up to 7mm in diameter

Electro Skin Cautery is used in Electrosurgery

Electrosurgery is done to destroy benign and malignant lesions, to control bleeding, and to cut or excise tissue.

Standard Vessel Sealing Cum Electrosurgical unit.

Graphical LCD Display.

400w Power.

99 programmable.

Patient Plate Conductivity monitoring system.

99 Programs Memory

Display : Graphical LCD Display

Dimension (W x H x D) : 360mm x 145 mm x 305 mm

Weight : 4.5 kg

Input Supply Voltage : 230 VAC + 15%@ 50/60Hz (110 V optional)

labycare

www.labycare.com | Email: labycare@gmail.com