

Labycare-Disposable Plates (SKSS-IE59)



Disposable Plates Manufacturer by LABYCARE

Description :

Material: 7 / 36 2core, ptfe insulated, silicon rubber insulated spc cable.

extra cable cord will be charged for this pad - silicon cable with 2jack pin one side with (t) junction or silicon cable with 2jack pin one the side with (t) junction and other side bipolar l&tconnector. (black color).

pad color: black

cable length: 3 meter.

jack pin size: 5 / 32 (4mm) or bipolar l&t connector for extra cable cord with this pad.

dimension: - bipolar l&t connector with cable cord with silicon pad.

design: - (universal designed).

packing: - sealed in high quality P. V. C. Plastic cover high quality branded look & safety.

Specifications :

Material: 7 / 36 2core, ptfe insulated, silicon rubber insulated spc cable.

extra cable cord will be charged for this pad - silicon cable with 2jack pin one side with (t) junction or silicon cable with 2jack pin one the side with (t) junction and other side bipolar l&tconnector. (black color).

pad color: black

cable length: 3 meter.

jack pin size: 5 / 32 (4mm) or bipolar l&t connector for extra cable cord with this pad.

dimension: - bipolar l&t connector with cable cord with silicon pad.

design: - (universal designed).

packing: - sealed in high quality P. V. C. Plastic cover high quality branded look & safety.

labycare

www.labycare.com | Email: labycare@gmail.com