

Labycare-Autoclave Pressure Cooker Powder Coated Electric (SKSSAS03)



Durable powder-coated electric autoclave ensures safe, efficient sterilization of medical instruments, glassware, and laboratory equipment. Easy-to-use, corrosion-resistant, and reliable, it delivers consistent high-pressure steam sterilization for hospitals and labs.

Description :

Model No: SKSSAS03

The Powder-Coated Electric Autoclave Pressure Cooker is an efficient sterilization device for laboratories, hospitals, and clinics. Its durable powder-coated body resists corrosion and ensures long-lasting use. Using high-pressure steam, it sterilizes instruments, glassware, and equipment safely and reliably. Equipped with temperature and pressure controls, it offers convenient, automated, and consistent sterilization for professional medical and laboratory applications.

Specifications :

Feature	Specification
Body Material	Powder-Coated Steel
Chamber Material	Stainless Steel
Capacity	10–50 liters
Heating Method	Electric
Sterilization Temperature	121°C – 134°C
Sterilization Pressure	15–30 psi
Control Panel	Digital / Manual
Safety Features	Pressure & Temperature Safety Valve, Overheat Protection
Operation	Fully Electric
Applications	Hospitals, Laboratories, Clinics

labycare

www.labycare.com | Email: labycare@gmail.com