

Labycare-Vertical Autoclave Single Wall Semi Automatic (SKSSAS013)



Vertical Autoclave Single Wall Semi Automatic Manufacturer by LABYCARE

Description :

Model No:- SKSSAS013

A vertical autoclave single wall semi-automatic model is a compact sterilizer widely used in clinics, laboratories, and research centers. Constructed with stainless steel, it features a single-walled chamber for efficient steam sterilization. The semi-automatic controls manage temperature and pressure, while safety valves and indicators ensure reliable operation and user protection.

Specifications :

Specification

Product Type

Material

Wall Type

Control System

Lid Locking

Temperature Range

Working Pressure

Capacity

Heating Source

Safety Features

Applications

Details

Vertical Autoclave Single Wall Semi Automatic

Stainless Steel (Corrosion Resistant)

Single Wall Construction

Semi-automatic temperature & pressure control

Heavy-duty locking mechanism

Up to 121°C

15–20 psi

5 to 30 liters (various models available)

Electric heating element

Pressure gauge, safety valve, auto cut-off

Hospitals, Clinics, Laboratories, Research Facilities

labycare

www.labycare.com | Email: labycare@gmail.com