

Labycare-Vertical Autoclave Single Wall Fully Automatic (SKSSAS014)



Vertical Autoclave Single Wall Fully Automatic Manufacturer by LABYCARE

Description :

Model No: SKSSAS014

The **Vertical Autoclave Single Wall Semi Automatic** is a robust and efficient sterilization unit designed for hospitals, laboratories, clinics, and research facilities. Built with high-quality stainless steel, it features a single-wall construction with a strong lid locking system for secure and airtight operation. The semi-automatic control ensures precise temperature and pressure regulation with minimal supervision, making sterilization safe and reliable. It is ideal for sterilizing surgical instruments, glassware, dressings, and laboratory materials. Compact and durable, this vertical autoclave offers uniform steam distribution, consistent performance, and ease of operation, making it a preferred choice in medical and laboratory environments.

Specifications :

Specification

Product Type

Material

Wall Type

Control System

Details

Vertical Autoclave Single Wall Semi Automatic

Stainless Steel (Corrosion Resistant)

Single Wall Construction

Semi-automatic temperature & pressure control

Lid Locking	Heavy-duty locking mechanism
Temperature Range	Up to 121°C
Working Pressure	15–20 psi
Capacity	5 to 30 liters (various models available)
Heating Source	Electric heating element
Safety Features	Pressure gauge, safety valve, auto cut-off
Applications	Hospitals, Clinics, Laboratories, Research Facilities

labycare

www.labycare.com | Email: labycare@gmail.com