

Labycare-Vertical Autoclave Double Wall Semi Automatic (SKSSAS017)



The Vertical Autoclave Double Wall Semi Automatic ensures safe, efficient sterilization with durable construction, semi-automatic controls, and reliable performance for hospitals and laboratories.

Description :

Model No : SKSSAS017

The **Vertical Autoclave Double Wall Semi Automatic** is a high-performance sterilizer designed for hospitals, laboratories, clinics, and pharmaceutical facilities. Built with durable stainless steel or mild steel double-wall construction, it provides excellent insulation, energy efficiency, and long-lasting performance. The semi-automatic control system ensures precise regulation of temperature and pressure with minimal supervision, delivering safe and consistent sterilization. Equipped with a heavy-duty lid locking system, safety valve, auto cut-off, and pressure gauge, it guarantees operator safety and reliability. Ideal for sterilizing surgical instruments, dressings, glassware, and laboratory materials, this autoclave combines durability, efficiency, and user-friendly operation for professional sterilization needs.

Specifications :

Specification

Product Type

Material

Wall Type

Control System

Heating Source

Temperature Range

Working Pressure

Locking System

Safety Features

Capacity

Applications

Details

Vertical Autoclave Double Wall Semi Automatic

Stainless Steel / Mild Steel

Double Wall (Improved insulation & durability)

Semi-automatic pressure & temperature control

Electric heating element

Up to 121°C

15–20 psi

Heavy-duty lid locking mechanism

Pressure gauge, safety valve, auto cut-off

5 to 100 liters (models available)

Hospitals, Clinics, Laboratories, Research & Pharma Facilities

labycare

www.labycare.com | Email: labycare@gmail.com