

Labycare-Vertical Autoclave Double Wall Fully Automatic (SKSS-AS018)



Vertical Autoclave Double Wall Fully Automatic Manufacturer by LABYCARE

Description :

Vertical Autoclave Double Wall Fully Automatic is made of Stainless Steel 304 Grade

Vertical Autoclave Full Automatic is top end model for sterilization of all critical application with exact reproducibility

It has Unique double wall design

Autoclave is fitted with Heavy duty jacket

Vertical Autoclave has fully automatic operation

Autoclave has μ based PID controller with dual display for Temperature & Time

Vertical Autoclave has reduced loading height for convenience.

The operating pressure is 1.2 / 2.1 kg / cm²

The operating temperature is 1210 C/1340 C

The Sterilization Time is 0-99 mins

The Voltage is 230V

The control of Vertical Autoclave is Digital

We are Manufacturer, Exporters and Suppliers of Vertical Autoclave Double Wall Fully Automatic at wholesale competitive prices from India.

Specifications :

Vertical Autoclave Double Wall Fully Automatic is made of Stainless Steel 304 Grade

Vertical Autoclave Full Automatic is top end model for sterilization of all critical application with exact reproducibility

It has Unique double wall design

Autoclave is fitted with Heavy duty jacket

Vertical Autoclave has fully automatic operation

Autoclave has μ p based PID controller with dual display for Temperature & Time

Vertical Autoclave has reduced loading height for convenience.

The operating pressure is 1.2 / 2.1 kg / cm²

The operating temperature is 1210 C/1340 C

The Sterilization Time is 0-99 mins

The Voltage is 230V

The control of Vertical Autoclave is Digital

We are Manufacturer, Exporters and Suppliers of Vertical Autoclave Double Wall Fully Automatic at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com