

Labycare-Horizontal Autoclave Cylindrical Deluxe (SKSSAS023)



Deluxe horizontal cylindrical autoclave ensures safe, uniform sterilization with stainless steel design, automatic control, and durable construction for hospitals.

Description :

Product Cat SKSSAS023

Specifications :

Parameter	Specification
Model Type	Horizontal Autoclave Cylindrical Deluxe
Chamber Construction	High-Grade Stainless Steel (SS 304 / SS 316 optional)
Outer Body	Polished Stainless Steel Double Wall
Chamber Type	Cylindrical Horizontal
Operation Type	Fully / Semi Automatic
Working Pressure	15–22 psi
Operating Temperature	Up to 134°C ±1°C
Heating System	Electric Heating / Steam Heating Option
Control Panel	Digital PID / Microprocessor Control with Timer
Door Type	Radial Locking with Safety Interlock
Safety Features	Automatic Pressure Cut-off, Low Water Cut-off, Safety Valve, Door Lock Mechanism
Pressure Gauge	0–60 psi Analog / Digital Display
Insulation	High Efficiency Glass Wool for Heat Retention
Finish	Deluxe Mirror Polish

Power Supply	440V AC, 3 Phase, 50 Hz
Mounting	Heavy-Duty Tubular Frame, Floor Mounted
Available Capacities	75 L, 150 L, 200 L, 300 L, 400 L, Custom Sizes
Applications	Hospitals, Pharma Industries, Research Labs, Bio Labs

labycare

www.labycare.com | Email: labycare@gmail.com