

Labycare-Table Top Autoclave Without Vacuum (SKSSAS028)



Compact stainless-steel tabletop autoclave ensures efficient, reliable sterilization without vacuum—ideal for dental, clinical, and laboratory applications.

Description :

Model No;:-SKSSAS028

The Table Top Autoclave Without Vacuum is a compact and efficient sterilization unit designed for clinics, laboratories, dental practices, and small hospitals. Constructed from high-quality stainless steel, it ensures durability and long-lasting performance. The autoclave operates using high-pressure saturated steam to sterilize instruments, glassware, and medical materials effectively. Its user-friendly design features a digital temperature controller, pressure gauge, safety valve, and timer for precise and reliable sterilization cycles. With its compact tabletop design and low maintenance requirements, this autoclave is ideal for small spaces requiring quick, efficient, and cost-effective sterilization.

Specifications :

Model Type	Table Top Autoclave Without Vacuum
Chamber Material	Stainless Steel (SS 304 / SS 316 optional)
Outer Body	Stainless Steel Polished Finish
Operation Type	Manual / Semi Automatic
Working Pressure	15–20 psi
Operating Temperature	Up to 121°C ±1°C
Sterilization Method	Steam Pressure Sterilization (Without Vacuum)
Heating System	Electric Heating Elements
Temperature Control	Digital PID Controller
Timer Range	0–60 Minutes
Pressure Gauge	Analog Display
Safety Features	Pressure Release Valve, Low Water Cut-off, Door

Power Supply
Capacity Options
Mounting Type
Finish
Applications

Locking System
220V AC, Single Phase, 50 Hz
10 L, 20 L, 25 L, 30 L, 40 L
Table Top
Mirror Polished Stainless Steel
Dental Clinics, Pathology Labs, Hospitals,
Research Institutes

labycare

www.labycare.com | Email: labycare@gmail.com