

Labycare-Table Top Autoclave With Vacuum (SKSS-AS029)



Table Top Autoclave With Vacuum Manufacturer by LABYCARE

Description :

Model No:- SKSSAS029

A **Table Top Autoclave with Vacuum** is a compact, high-efficiency sterilizer designed for medical, dental, and laboratory use. It uses pre and post vacuum cycles to ensure complete air removal and thorough steam penetration, providing fast, reliable sterilization and effective drying for hollow, wrapped, or porous instruments.

Specifications :

Parameter	Details
Capacity	23 L (247 mm x 470 mm)
Sterilization Class	Class B (EN 13060)
Temperature Options	121 °C and 134 °C
Sterilization Cycle Time	Approx. 20 – 25 min
Maximum Pressure	~2.3 bar
Maximum Vacuum Depth	~0.9 bar (negative)
Heating	Built-in Steam Generator
Drying	Post-vacuum drying system
Programs	Multiple programs (wrapped, unwrapped, B&D test, vacuum test, etc.)
Control System	Digital microprocessor with PLC logic
Display & Logging	LCD with USB data logging
Safety Features	Automatic pressure/temperature cutoff, door interlock, over-pressure valve

Water Tank Capacity	~3.5 L
Water Consumption per Cycle	0.16 – 0.18 L
Power Supply	110 V / 60 Hz
Electrical Power	~1950 W
Exterior Dimensions	495 x 700 x 410 mm (W x D x H)
Material	SUS 304 Stainless Steel Chamber
Additional Features	B&D & Helix tests, Printer/USB port, Safety alarms

labycare

www.labycare.com | Email: labycare@gmail.com