

Labycare- Dental Autoclave with Accessories (SKSSVA355G)



Dental Autoclave with Accessories Manufacturer by LABYCARE

Description :

Model No:- SKSSVA355G

A **Dental Autoclave** is an essential piece of equipment used in dental clinics and laboratories for sterilizing instruments, ensuring that all dental tools and devices are free from bacteria, viruses, and other pathogens. The autoclave uses high-pressure steam to achieve sterilization, providing an efficient and reliable method for infection control. This unit is designed with user-friendly controls and is equipped with essential safety features, making it suitable for high-volume dental practices. The **Dental Autoclave with Accessories** comes complete with the necessary tools for effective operation, including trays, baskets, and tongs, for safely placing and handling instruments during the sterilization process.

Specifications :

Parameter	Details
Product Type	Dental Steam Autoclave with Accessories
Capacity	18 L / 23 L
Material	Stainless Steel (SUS 304)
Sterilization Temperature	121°C / 134°C
Pressure Range	1.1 – 2.3 bar
Cycle Duration	20 – 30 minutes
Operation Type	Semi / Fully Automatic
Drying System	Automatic or Vacuum Drying
Included Accessories	Instrument trays, holders, tubing, water reservoir, sealing gasket

Display Type
Power Supply
Safety Features

Applications

Digital / LCD Panel
220 V AC, 50 Hz
Door interlock, over-temperature & over-pressure protection
Dental clinics, orthodontic centers, oral surgery, hospitals

labycare

www.labycare.com | Email: labycare@gmail.com