

Labycare-Table Top Front Loading Autoclave (SKSS1101A)



Table Top Front Loading Autoclave Manufacturer by LABYCARE

Description :

Product Cat # SKSS1101A

The Table Top Front Loading Autoclave is a compact, high-efficiency steam sterilizer designed for medical, dental, and laboratory use. Its front-loading design ensures easy access, while digital controls and automatic cycles provide fast, reliable sterilization and drying of instruments, maintaining the highest hygiene and safety standards in clinics and labs.

Specifications :

Feature	Specification
Model Type	Table Top Front Loading Autoclave
Chamber Volume	18L / 23L / 25L
Chamber Material	SUS 304 / 316 Stainless Steel
Outer Body	Powder-Coated / Stainless Steel
Sterilization Temp	121°C & 134°C
Max Pressure	2.3?bar
Vacuum System	Optional Pre & Post Vacuum
Cycle Time	20–30?minutes
Control Panel	Digital / LCD with Programmable Cycles
Drying	Automatic / Vacuum Drying
Water Tank	3.5?L Built-In

Power
Safety Features
Applications
Optional Accessories

220?V AC, 50?Hz, 1950?W
Door Interlock, Over-Temp & Pressure Protection
Hospitals, Dental Clinics, Labs, Veterinary Centers
Instrument Trays, Holders, USB/Printer Interface

labycare

www.labycare.com | Email: labycare@gmail.com