

Labycare-Bowl and Utensil Sterilizer (SKSSAS033)



Bowl and Utensil Sterilizer Manufacturer by LABYCARE

Description :

Model No:- SKSSAS033

The **Bowl and Utensil Sterilizer** is a stainless-steel electric unit designed to sterilize medical bowls, trays, and instruments using steam. It ensures complete disinfection, featuring thermostatic control, water outlet, and durable construction. Ideal for hospitals, clinics, and laboratories to maintain high hygiene standards and prevent cross-contamination.

Specifications :

**Specification
Material**

Details

Stainless Steel (SS 304) – corrosion-resistant and

Sterilization Method	durable
Operation	Steam sterilization using electric heating element
Temperature Range	Electric with thermostatic temperature control
Power Supply	Ambient +5°C to 99°C
Control Panel	220V / 50Hz AC single phase
Safety Features	Digital timer and temperature controller
Lid Mechanism	Low-water level cut-off and automatic shut-off
Drainage System	Pedal-operated lid for hands-free use
Tray Stand	Built-in drain outlet for easy water removal
Available Sizes	Optional stainless steel tray stand with hanger
	20"×16"×16", 20"×20"×20", 24"×20"×20", 24"×24"×20"
Applications	Used for sterilizing bowls, trays, basins, and surgical instruments in hospitals and clinics

labycare

www.labycare.com | Email: labycare@gmail.com