

Labycare-Infant Radiant Warmer Silver without Phototherapy (SKSS-BCE02)



Infant Radiant Warmer Silver without Phototherapy Manufacturer by Labycare

Description :

The Infant Radiant Warmer Silver without phototherapy is designed to maintain a stable thermal environment for newborns. It delivers uniform radiant heat with precise temperature control, ensuring safe, comfortable warming and unobstructed access for clinical procedures in NICUs and delivery rooms.

Specifications :

Type: Infant Radiant Warmer (without phototherapy)

Finish: Silver

Heat Source: Overhead radiant heater

Temperature Control: Servo and manual modes

Display: Digital temperature indicator

Mattress: Soft, X-ray permeable mattress

Alarms: Audio-visual alarms for safety

Power Supply: 220–240 V AC, 50 Hz

Application: NICU, labor room, neonatal care units

labycare

www.labycare.com | Email: labycare@gmail.com