

Labycare-Infant Radiant Warmer Gold (SKSS-BCE03)



Infant Radiant Warmer Gold Manufacturer by LABYCARE

Description : The Infant Radiant Warmer Gold provides advanced thermal care for newborns with reliable radiant heating and precise temperature control. Designed for intensive neonatal care, it ensures uniform warmth, quick access to the infant, enhanced safety features, and a premium gold-finish structure suitable for modern NICUs.

Specifications :

Type: Infant Radiant Warmer

Finish: Gold

Heat Source: Overhead radiant heater

Temperature Control: Servo and manual modes

Display: Digital temperature display

Mattress: X-ray permeable, washable mattress

Alarms: Audio-visual safety alarms

Power Supply: 220–240 V AC, 50 Hz

Application: NICU, delivery rooms, neonatal care units

