

Labycare-Infant Radiant Warmer Gold without Phototherapy (SKSS-BCE04)



Infant Radiant Warmer Gold without Phototherapy Manufacturer by LABYCARE

Description : The Infant Radiant Warmer Gold without phototherapy is designed for effective thermal regulation of newborns. It provides consistent radiant heat, precise temperature control, and easy infant access. The premium gold finish, safety alarms, and ergonomic design make it ideal for NICUs and delivery rooms

Specifications :

Type: Infant Radiant Warmer (without phototherapy)

Finish: Gold

Heat Source: Overhead radiant heater

Temperature Control: Servo and manual modes

Display: Digital temperature indicator

Mattress: X-ray permeable, washable

Alarms: Audio-visual alarms for temperature and power failure

Power Supply: 220–240 V AC, 50 Hz

Application: NICU, labor room, neonatal care units

