

Labycare-Infant Radiant Warmer Diamond (SKSS-BCE05)



Infant Radiant Warmer Diamond Manufacturer by LABYCARE

Description :

The Infant Radiant Warmer Diamond is an advanced neonatal care device designed to provide precise thermal regulation for newborns. With uniform radiant heating, accurate temperature control, and enhanced safety systems, it ensures optimal warmth, comfort, and easy clinical access in NICUs and delivery rooms.

Specifications :

Type: Infant Radiant Warmer

Model: Diamond

Heat Source: Overhead radiant heater

Temperature Control: Servo and manual modes

Display: Digital temperature display

Mattress: X-ray permeable, washable mattress

Alarms: Audio-visual alarms for safety and faults

Power Supply: 220–240 V AC, 50 Hz

Application: NICU, delivery rooms, neonatal care units

labycare

www.labycare.com | Email: labycare@gmail.com