

Labycare-Infant Radiant Warmer Diamond without Phototherapy (SKSS-BCE06)



Infant Radiant Warmer Diamond without Phototherapy Manufacturer by LABYCARE

Description :

The Infant Radiant Warmer Diamond without phototherapy is designed for precise thermal care of newborns. It delivers uniform radiant heat with accurate temperature control, ensuring infant safety and comfort. The Diamond model offers enhanced reliability, clear digital monitoring, and easy access for clinical procedures in NICUs.

Specifications :

Type: Infant Radiant Warmer (without phototherapy)

Model: Diamond

Heat Source: Overhead radiant heater

Temperature Control: Servo and manual modes

Display: Digital temperature indicator

Mattress: X-ray permeable, washable

Alarms: Audio-visual safety and power failure alarms

Power Supply: 220–240 V AC, 50 Hz

Application: NICU, delivery rooms, neonatal care units

labycare

www.labycare.com | Email: labycare@gmail.com