

Labycare-Infant Radiant Warmer Diamond with Drawer & Phototherapy (SKSS-BCE08)



Infant Radiant Warmer Diamond with Drawer & Phototherapy Manufacturer by LABYCARE

Description :

The Infant Radiant Warmer Diamond with drawer and phototherapy offers comprehensive neonatal care. It combines uniform radiant heating with effective LED phototherapy for jaundice management. The integrated storage drawer improves workflow, while precise controls and safety alarms ensure reliable, comfortable, and safe care in NICUs and delivery rooms.

Specifications :

- Type: Infant Radiant Warmer with Phototherapy
- Model: Diamond
- Features: Built-in drawer, integrated phototherapy unit
- Heat Source: Overhead radiant heater
- Phototherapy: High-intensity LED blue light
- Temperature Control: Servo and manual modes

- Display: Digital temperature display
- Mattress: X-ray permeable, washable
- Alarms: Audio-visual alarms for safety and faults
- Power Supply: 220–240 V AC, 50 Hz
- Application: NICU, delivery rooms, neonatal care units

labycare

www.labycare.com | Email: labycare@gmail.com