

## Labycare-Drug Trolley (SKSS835A)



Drug Trolley

### **Description :**

Drug Trolley

### **Specifications :**

#### **Drug Trolley – Technical Specifications**

- **Material:** Stainless steel (SS304) / Epoxy-coated mild steel
- **Dimensions:** 700 x 450 x 900 mm
- **Drawers:** 2-4 lockable, smooth-gliding drawers
- **Shelves:** 1-2 with raised guardrails
- **Casters:** 4 swivel casters (100 mm diameter, 2 with brakes)
- **Handles:** Ergonomic push handle
- **Load Capacity:** Up to 100 kg
- **Finish:** Corrosion-resistant, easy-to-clean
- **Optional Accessories:** IV pole, waste bin, writing board
- **Locking:** Centralized locking system
- **Warranty:** 1-year against manufacturing defects

---

labycare

www.labycare.com | Email: labycare@gmail.com