

Labycare-Open Care System Silver Without Phototherapy (SKSS-BCE10)



Open Care System Silver Without Phototherapy

Description :

The Open Care System Silver without phototherapy is a neonatal radiant warmer designed to maintain optimal body temperature for newborns. It provides uniform heating, unobstructed access for medical care, and precise temperature control, making it suitable for NICUs, delivery rooms, and neonatal care units.

Specifications :

Type: Neonatal open care radiant warmer

Phototherapy: Not included

Heater: Overhead radiant heating element

Temperature control: Servo-controlled and manual modes

Display: Digital temperature indicator with alarms

Baby tray: X-ray translucent mattress

Alarms: Audio-visual alarms for temperature deviation and power failure

Power supply: 220–240 V AC, 50 Hz
