

Labycare-Open Care System Gold (SKSS-BCE011)



Open Care System Gold Manufacturer by LABYCARE

Description :

Open Care System Gold is a neonatal open care radiant warmer that provides advanced thermal support for newborns. It features a microprocessor-controlled heating system with LCD display, audible/visual alarms, and a spacious bassinet for easy access. The heater can swivel for procedures and ensures continuous operation even on probe failure.

Specifications :

- Control Modes: Manual, Servo (baby) control with automatic conversion on probe failure.
- Display: Large LCD showing skin/air and set temperatures; °C/°F selectable.
- Heater: 650 W silica/quartz/ceramic radiant heater with integral examination light; swivel-able for X-rays/procedures.
- Alarms: High/Low temp, system fail, probe failure, power fail (audible & visual).
- Baby Bed: Spacious with drop-down/fold-down sides for access.
- Mobility: Castors for transport
- Power: ~200 V/50 Hz, fuse protection

labycare

www.labycare.com | Email: labycare@gmail.com