

Labycare-Open Care System Diamond (SKSS-BCE013)



Open Care System Diamond Manufacturer by LABYCARE

Description :

The **Open Care System Diamond** is a neonatal care unit providing precise temperature control through automatic, manual, and servo modes. It features a large LCD, safety alarms, and integrated LED phototherapy. Designed for effective infant warming, monitoring, and comfort in medical settings, it enhances care efficiency and safety.

Specifications :

- Dimensions: 185 cm (H) × 105 cm (L) × 65 cm (W)
- Power: 200 V, 50 Hz, 650 W
- Control Modes: Automatic, Manual, Servo
- Display: Large LCD
- Safety Alarms: High/Low temperature, system alerts
- Extras: LED phototherapy, timer, pre-warm function

labycare

www.labycare.com | Email: labycare@gmail.com