

Labycare-Electric Sterilizer (SKSSES908)



Electric Sterilizer

Description :

Electric Sterilizer

Specifications :	Electric Sterilizer	Electric Sterilizer 12 x 6 x 5.5 inch.
	Brand	SKSS
	- Instrument Sterilizers, electric, with hinged lid fitted with immersion rod and water outlet tap - Seamless stainless steel construction. - Fitted with automatic cutoff safety pin device, thermostat &	

tray lifting
arrangement.

labycare

www.labycare.com | Email: labycare@gmail.com