

## Labycare-Baby Infant Incubator With Phototherapy (SKSS-BCE022)



Baby Infant Incubator With Phototherapy

### **Description :**

A Baby Infant Incubator with Phototherapy combines a controlled thermal environment with integrated blue-light therapy to treat neonatal jaundice. It ensures stable temperature, humidity, and oxygen support while delivering safe, effective phototherapy for premature and newborn infants in NICUs.

### **Specifications :**

Temperature control range: 20°C–37°C

Humidity control: 40%–90% RH

Integrated LED phototherapy unit (blue light 430–490 nm)

Microprocessor-based control panel

Transparent acrylic hood

Adjustable mattress with X-ray tray

Audible and visual safety alarms

Power supply: 220–240 V AC

