

Labycare-Penguin Sucker (SKSS-BCE033)



Penguin Sucker Manufacturer by LABYCARE

Description :

Penguin Sucker is a neonatal mucus suction device used to gently clear secretions from a newborn's mouth and nose. Designed for safe, effective aspiration, it ensures clear airways immediately after birth and during routine neonatal care.

Specifications :

- Soft, medical-grade, non-toxic PVC
- Gentle suction with controlled pressure
- Clear chamber for visual monitoring
- Rounded tip to prevent tissue injury
- Easy to use, clean, and sterilize
- Suitable for neonatal and infant use

labycare

www.labycare.com | Email: labycare@gmail.com