

## Labycare-Fetoscope Plastic (SKSS-BCE038)



Fetoscope Plastic

### **Description :**

A Plastic Fetoscope is a simple, durable instrument used by healthcare professionals to listen to fetal heart sounds during pregnancy. Its non-electric, acoustic design ensures reliable performance, easy handling, and suitability for routine prenatal checkups in hospitals, clinics, and outreach programs.

### **Specifications :**

Material: High-quality medical-grade plastic

Type: Acoustic fetoscope

Use: Fetal heart sound monitoring

Design: Pinard type / trumpet shape

Length: Approx. 18–20 cm

Operation: Non-electric, manual

Application: Antenatal care, obstetrics

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)