

Labycare-Fetoscope Wooden (SKSS-BCE039)



Fetoscope Wooden Manufacturer by LABYCARE

Description :

A Wooden Fetoscope is a traditional acoustic instrument used to monitor fetal heart sounds during pregnancy. Made from smooth, polished wood, it provides clear sound transmission without electricity, making it reliable for routine antenatal examinations in hospitals, clinics, and community healthcare settings.

Specifications :

Material: Polished natural wood

Type: Acoustic fetoscope (Pinard type)

Use: Fetal heart sound monitoring

Length: Approx. 18–20 cm

Operation: Non-electric, manual

Finish: Smooth, splinter-free surface

Application: Obstetrics and antenatal care

labycare

www.labycare.com | Email: labycare@gmail.com