

Labycare-Fetoscope Aluminium (SKSS-BCE40)



Fetoscope Aluminium Manufacturer by LABYCARE

Description :

An Aluminum Fetoscope is a durable acoustic instrument used to listen to fetal heart sounds during pregnancy. Its lightweight metal construction ensures excellent sound conduction and long service life, making it suitable for routine antenatal examinations in hospitals, clinics, and maternity care centers.

Specifications :

Material: Medical-grade aluminium

Type: Acoustic fetoscope (Pinard type)

Use: Fetal heart auscultation

Length: Approx. 18–20 cm

Operation: Non-electric, manual

Finish: Smooth, corrosion-resistant

Application: Obstetrics and antenatal care

labycare

www.labycare.com | Email: labycare@gmail.com