

Labycare-Infant Vein Viewer (SKSS-BCE042)



Infant Vein Viewer

Description :

An Infant Vein Viewer is a non-invasive medical device that helps clinicians clearly visualize veins in infants and neonates. It uses safe infrared light to project vein patterns onto the skin, improving accuracy during venipuncture, reducing pain, and minimizing repeated needle attempts.

Specifications :

Uses near-infrared (NIR) light technology

Safe, non-invasive, and radiation-free

Real-time vein projection on skin

Adjustable brightness and focus

Lightweight, portable, handheld design

Suitable for infants and neonates

labycare

www.labycare.com | Email: labycare@gmail.com