

Labycare-Infrared Vein Finder Multi Purpose (SKSS-BCE044)



Infrared Vein Finder Multi Purpose

Description :

An Infrared Vein Finder is a non-invasive medical device that uses near-infrared technology to detect and display veins beneath the skin. It assists healthcare professionals in accurate venipuncture across all age groups, reducing procedure time, patient discomfort, and repeated needle attempts.

Specifications :

Near-infrared (NIR) vein detection technology

Clear real-time vein visualization

Non-invasive, radiation-free, and safe

Adjustable brightness and imaging modes

Portable, handheld, or stand-mounted options

Suitable for pediatric, adult, and geriatric use

labycare

www.labycare.com | Email: labycare@gmail.com