

Labycare-Oxygen Hood With Nozzle (SKSS-BCE046)



Oxygen Hood With Nozzle Manufacturer by LABYCARE

Description :

An Oxygen Hood with Nozzle is a medical oxygen delivery device used primarily for infants and neonates. It provides a controlled, high-concentration oxygen environment while ensuring patient comfort. The integrated nozzle allows easy connection to an oxygen source for safe and efficient respiratory support.

Specifications :

Made of transparent, medical-grade plastic

Designed for neonatal and pediatric use

High oxygen concentration delivery capability

Smooth edges for patient comfort and safety

Integrated nozzle for oxygen inlet connection

Easy to clean and reusable or disposable options

labycare

www.labycare.com | Email: labycare@gmail.com