

Labycare-Fetal Heart Doppler (SKSS-BCE048)



Fetal Heart Doppler Manufacturer by LABYCARE

Description :

A Fetal Heart Doppler is a portable ultrasound device used to detect and monitor fetal heart rate during pregnancy. It provides clear audio and visual indications of fetal heartbeat, assisting healthcare professionals in routine antenatal examinations and ensuring fetal well-being in a safe, non-invasive manner.

Specifications :

High-sensitivity ultrasound Doppler technology

Accurate fetal heart rate measurement

Clear audio output with built-in speaker

Digital display for heart rate reading

Safe, non-invasive, and easy to use

Suitable for clinical and prenatal care use

