

Labycare-Autopsy Table Stainless Steel (SKSS378A)



Autopsy Table Stainless Steel Manufacturer by LABYCARE

Description :

Stainless steel autopsy table is designed for hygienic and efficient post-mortem examinations. It features corrosion-resistant SS construction, smooth seamless surface, integrated drainage system, raised edges, and sturdy frame. Easy to clean, disinfect, and maintain, ensuring safety, durability, and compliance with mortuary standards.

Specifications :

Material: High-grade stainless steel (SS 304 / SS 316 optional)

Construction: Seamless, corrosion-resistant, hygienic finish

Top: Sloped surface with raised edges

Drainage: Integrated drain outlet with strainer

Frame: Heavy-duty tubular stainless steel

Load Capacity: Suitable for adult cadavers

Finish: Smooth, matte or mirror finish

Cleaning: Easy to wash and disinfect

Application: Hospitals, mortuaries, forensic labs

labycare

www.labycare.com | Email: labycare@gmail.com