

Labycare-Autopsy Table Electro Hydraulic (SKSS961)



Autopsy Table Electro Hydraulic Manufacturer by LABYCARE

Description : Electro hydraulic autopsy table is designed for advanced post-mortem procedures with smooth height adjustment and tilting. It ensures ergonomic working conditions, superior hygiene, and safe handling of cadavers. Ideal for modern mortuaries and forensic departments requiring precision, efficiency, and operator comfort.

Specifications :

Operation: Electro-hydraulic, push-button control

Material: Stainless steel SS 304 / SS 316

Height Adjustment: Smooth vertical movement

Table Top: Seamless, sloped with raised edges

Drainage: Built-in drain with strainer

Tilting: Head-up / foot-up tilt facility

Load Capacity: Suitable for adult bodies

Power Supply: Standard AC electric supply

Finish: Corrosion-resistant, easy-clean surface

Application: Hospitals, mortuaries, forensic labs
