

Labycare-Autopsy Forensic Workstation (SKSS-AP03)



Autopsy Forensic Workstation Manufacturer by LABYCARE

Description : An autopsy forensic workstation is a specialized unit designed for safe, hygienic, and efficient post-mortem examination. Built with high-grade stainless steel, it integrates ventilation, drainage, and lighting systems. The workstation ensures contamination control, ergonomic workflow, and smooth instrument handling, supporting detailed forensic analysis in hospitals, labs, and mortuaries.

Specifications :

Material: Stainless Steel SS 304 heavy-duty construction

Dimensions: Approx. 2100–2400 mm (L) × 800–900 mm (W) × 900–1000 mm (H)

Work Surface: Sink, perforated top, integrated drainage

Ventilation: Built-in downdraft/side draft exhaust system

Utilities: Hot & cold water, spray gun, faucet

Lighting: Overhead LED examination lighting

Waste Management: Removable waste collector, drain traps

Electrical: Waterproof sockets, control panel

Additional Features: Instrument shelf, splash guard, odor control system

Optional: Height adjustment, camera mount, refrigeration connection

labycare

www.labycare.com | Email: labycare@gmail.com