

## Labycare-Dead Body Transfer Trolley (SKSSMC2D)



### Dead Body Transfer Trolley

**Description :** A dead body transfer trolley is designed for safe, smooth, and dignified transportation of bodies in hospitals, mortuaries, and emergency facilities. Built with a sturdy stainless-steel frame, it features smooth-rolling wheels, secure side rails, and an easy-to-clean surface, ensuring durability, hygiene, and reliable handling during body transfer operations

### Specifications :

Frame Material: Stainless Steel (SS 304/202) or powder-coated steel

Tray/Top: Removable or fixed stainless-steel stretcher

Dimensions: Approx. 1900–2000 mm (L) × 550–650 mm (W) × 800–900 mm (H)

Load Capacity: 150–250 kg

Wheels: 4 heavy-duty castor wheels (2 with brakes)

Handles: Push/pull handles for smooth movement

Finish: Corrosion-resistant, easy-to-clean polished surface

Construction: Welded, durable, and stable frame

Additional Options: Collapsible design, draped cover, locking mechanism, central braking system

