

Labycare-Dead Body Wash Trolley (SKSSMC2C)



Dead Body Wash Trolley Manufacturer by LABYCARE

Description :

A dead body wash trolley is designed for hygienic and efficient washing of bodies in hospitals and mortuaries. Made with a durable stainless-steel frame, it includes a perforated top or drain system for easy water flow. Heavy-duty wheels, adjustable height options, and easy-clean surfaces ensure safe and convenient handling.

Specifications :

Material: Stainless Steel SS 304/202 construction

Top Surface: Perforated or slotted SS top for smooth drainage

Dimensions: Approx. 1900–2000 mm (L) × 600–700 mm (W) × 800–900 mm (H)

Load Capacity: 150–250 kg

Drainage: Built-in drain outlet with slope design

Frame: Rigid welded frame with corrosion resistance

Wheels: 4 heavy-duty castors (2 with brakes)

Finish: Polished, easy-to-clean surface

Optional Features: Adjustable height, detachable top, water collection tray, side rails

labycare

www.labycare.com | Email: labycare@gmail.com