

Labycare-Mortuary Trolley Height Adjustable (SKSSMC2A)



Mortuary Trolley Height Adjustable Manufacturer by LABYCARE

Description : A height-adjustable mortuary trolley is designed for efficient body handling in hospitals and mortuaries. Its adjustable mechanism allows smooth elevation for examination, transfer, or storage. Built with a strong stainless-steel frame, it ensures hygiene, stability, and durability. Heavy-duty wheels and an easy-to-clean surface provide safe, controlled, and convenient operation.

Specifications :

Material: Stainless Steel SS 304/202

Height Adjustment: Manual or hydraulic; approx. 600–900 mm

Dimensions: 1900–2000 mm (L) × 600–700 mm (W)

Top: Removable or fixed SS stretcher/slab

Load Capacity: 150–250 kg

Wheels: 4 heavy-duty castors (2 with brakes)

Frame: Rigid, corrosion-resistant, welded construction

Additional Options: Central braking system, locking mechanism, drainage tray, side rails

labycare

www.labycare.com | Email: labycare@gmail.com