

## Labycare-Dead Body Concealment Trolley (SKSSMC2B)



Dead Body Concealment Trolley Manufacturer by LABYCARE

**Description :** A dead body concealment trolley is designed for discreet and dignified transport of bodies within hospitals and mortuaries. It features a covered structure with a stainless-steel stretcher, smooth-rolling wheels, and a durable frame. The concealed design maintains privacy, while its easy-to-clean surfaces ensure hygiene and efficient movement.

### **Specifications :**

Material: Stainless Steel SS 304/202 frame

Cover: Fabric/PU-coated/zippered concealment canopy

Dimensions: Approx. 1900–2000 mm (L) × 600–700 mm (W) × 900–1100 mm (H)

Top: Removable stainless-steel stretcher/slab

Load Capacity: 150–250 kg

Wheels: 4 heavy-duty castors (2 with brakes)

Frame: Fully welded, corrosion-resistant construction

Features: Concealment curtain/canopy, hygienic surface, push handles

Optional: Central braking, adjustable height, waterproof cover

