

Labycare-Organ Weighing Scale (SKSS-AP014)



Organ Weighing Scale Manufacturer by LABYCARE

Description : Organ weighing scale is designed for precise measurement of human organs during autopsy and pathology procedures. It offers high accuracy, hygienic stainless steel construction, and easy-to-read digital display. Compact, reliable, and easy to clean, it ensures accurate documentation in medical, forensic, and laboratory environments.

Specifications :

Type: Digital / Electronic

Capacity: 5–30 kg (as required)

Accuracy: High-precision sensor

Platform: Stainless steel, corrosion resistant

Display: Clear LCD/LED digital display

Functions: Tare, zero, auto calibration

Power Supply: Battery operated / AC adapter

Design: Compact, stable, easy to clean

Application: Autopsy rooms, pathology labs, hospitals
