

Labycare-Embalming Workstation (SKSS-AP015)



Embalming Workstation Manufacturer by LABYCARE

Description :

An Embalming Workstation is a complete, hygienic setup designed for professional embalming procedures. Constructed from corrosion-resistant stainless steel, it ensures durability, easy cleaning, and efficient workflow with integrated drainage, storage, and fluid management features.

Specifications :

Material: SS 304 stainless steel

Worktop: Sloped with drainage outlet

Sink: Integrated with water inlet

Storage: Cabinets / drawers

Plumbing: Hot & cold water provision

Finish: Matte / mirror polish

Application: Mortuaries, hospitals, anatomy labs

labycare

www.labycare.com | Email: labycare@gmail.com