

Labycare-Oscillating Autopsy Saw (SKSS-AP016)



Oscillating Autopsy Saw Manufacturer by LABYCARE

Description : An Oscillating Autopsy Saw is designed for precise, controlled cutting of bone during post-mortem examinations. Its oscillating action minimizes tissue damage and aerosol generation, ensuring safety, accuracy, and efficiency for forensic and pathological procedures.

Specifications :

Type: Electric / Pneumatic oscillating saw

Blade: Stainless steel, interchangeable

Speed: Variable oscillation control

Handle: Ergonomic, non-slip grip

Safety: Low vibration, reduced splatter

Power Supply: 220–240 V AC (electric)

Application: Autopsy rooms, forensic labs

labycare

www.labycare.com | Email: labycare@gmail.com